



MetScape 2.3.2

User Manual

A Plugin for Cytoscape

National Center for Integrative Biomedical Informatics

July 2012

Contents

OVERVIEW	4
About data sources	4
Workflow overview	4
INSTALLING CYTOSCAPE AND THE METSCAPE PLUGIN	5
Launching Cytoscape and the MetScape plugin from the Web	5
Installing Cytoscape and the MetScape plugin on your local computer	5
ENTERING DATA	6
Option 1: Enter a list of compounds	7
Option 2: Load an experiment file	9
VISUALIZATION	14
METSCAPE TAB	15
RULES USED TO BUILD DIFFERENT NETWORK TYPES	16
DATA PANEL DISPLAY	17
Choosing attributes provided through MetScape	17
Rearranging attributes in the data panel	20
Pathway filter in the data panel	21
Concept filter in the data panel	22
METSCAPE RESULTS PANEL (ADDITIONAL NODE INFORMATION)	25
ANIMATION OF DATA	26
Building an Animation	26
Manipulating the animation: zooming	28
Manipulating the animation: bar chart	29
Playing the animation	30
MANIPULATING NETWORKS	30
Expanding a network	30
Option 1: Expanding a network in the current network window (Example using Compound- Reaction Network)	30
Option 2: Expanding a network in a new window (Example using Compound Network from Selected Genes)	32
Collapsing a network	33
Creating a subnetwork	33
Destroying a network	35

SAVING AND REOPENING A SESSION 35
 Saving a session 35
 Reopening a session 35

Please note that due to continuous software upgrades, the images in this handout may not exactly mimic what you see on the screen.

OVERVIEW

About data sources

MetScape is a plug-in for [Cytoscape](#), the bioinformatics network visualization tool. The plug-in can be used to visualize and interpret metabolomics and gene expression data in the context of human metabolic networks.

MetScape uses a metabolite database developed by extracting and integrating information from the following sources:

1. Edinburgh Human Metabolic Network (EHMN) —
<http://www.ehmn.bioinformatics.ed.ac.uk/>
2. KEGG COMPOUND Database —
<http://www.genome.jp/kegg/compound/>

MetScape allows users to load a list of metabolites with experimentally determined concentrations, a list of genes with experimentally determined expression values, and a list of concepts or pathways, and display them in the context of relevant metabolic networks.

Workflow overview

With MetScape, you can:

- Trace the connections between metabolites and genes.
- Integrate multidimensional data
- Visualize compound, reaction, enzyme and gene networks and display compound structures, as well as information for reactions, enzymes, genes, and pathways.
- Visually animate changes in compound concentrations over time and across experimental conditions.

The basic steps in the workflow include:

1. [Enter data](#). You can type or paste a list of compounds and/or genes, load a file containing experimental data, or start from a biological pathway.
2. [Select compound and reaction attributes](#). Choose which attributes to display in a table as you work with your visual network graph.
3. Explore the Visual network and table of attributes:
 - [Expand](#) and [collapse](#) a network.
 - [Create a subnetwork](#)
 - Visualize your data in a wide variety of network layouts provided by Cytoscape.
 - Use color, size, and other effects to visually reflect a set of attribute values.
4. [Save your session and reopen it later](#).

INSTALLING CYTOSCAPE AND THE METSCAPE PLUGIN

You can use the MetScape plugin for Cytoscape in one of two ways:

- **By launching Cytoscape and the MetScape plugin from the [MetScape Webpage](#).** With this method, you launch a temporary copy of Cytoscape and the MetScape plugin. Use this method when you don't have the administrative privileges to install software on your computer. You will not have access to other plugins when you use this method.
- **By installing Cytoscape and the MetScape plugin on your local computer.** With this method, you can use other plugins that work with MetScape data, such as Shortest Path and MCODE.

Launching Cytoscape and the MetScape plugin from the Web


To launch Cytoscape and the MetScape plugin from the MetScape Webpage:


1. Go to <http://metscape.ncibi.org/metscape2/>.
2. On the **Webpage**, click the **Launch Cytoscape w/ MetScape 2 via Java Web Start** link.

After launching the Web version, you can do anything described in this User Guide.

You can save your Cytoscape session on your local computer and reopen it at a later time after launching the application from the MetScape Webpage again.

Installing Cytoscape and the MetScape plugin on your local computer

1. Install Cytoscape on your computer. For more information, go to the Cytoscape website at <http://cytoscape.org/>.
2. After Cytoscape is installed, start the application.
3. To install the MetScape plugin, select **Plugins -> Manage Plugins** from the Cytoscape menu.
4. Under **Available to Install**, click  (the plus sign) to expand the **Integrated Analysis** group.



5. Find and click the most recent version of the MetScape plugin.
6. Click **Install**.
7. When installation is complete, click **Close**.

A MetScape option is added to the Cytoscape **Plugins** menu.

Registration: MetScape is a free program. We ask you to register because it helps us to keep track of the number of downloads. Your information will be stored in a secure database and we will not share it with anyone. We may send you infrequent e-mails about future MetScape releases. You may also choose to decline the registration. If you decide to register later, the registration window can be found on the main MetScape menu.


Note: This plugin requires Cytoscape 2.8 to run correctly.

ENTERING DATA


1. To begin a Cytoscape session with the MetScape plugin, first start Cytoscape.
2. Choose one of the following methods to get started.

Option 1: Enter a list of compounds


1. Select **Plugins -> MetScape-> Build Network** from the Cytoscape menu.
2. A MetScape tab now appears on the left of the Cytoscape screen.
3. Select a species. Choices are: human, rat, mouse
4. Manually enter (or copy and paste) compound ID(s) or name(s) and/or Entrez gene ID(s) or symbol(s). Click **Add** and enter the appropriate ID type in the popup box. Lists of IDs should be separated by spaces.



5. Click **OK**
6. If you enter a compound name, it will map to its KEGG id; a popup window will appear. If there is more than one potential match, use the dropdown arrow to choose the best match. If the compound is not found in the system, it will say "Not Found."




- Click **OK**
- After the experimental data has been loaded, if any of the genes/compounds that you submitted were not mapped to the database objects, a **Missing Data** window will appear. After viewing the **Missing Data** information, click **SAVE** to save the data for later viewing or click **OK** to close the window without making it available for future viewing. (See the Missing Data Window section below, under Option2: Load an experiment file, for more information about missing data).




- Select a **Network Type** by selecting one of the following from the dropdown menu:


- **Compound-Reaction-Enzyme-Gene** (1)
- **Compound-Reaction** (2)
- **Compound-Gene** (3)
- **Compound** (4)




1



2



3



4

Note: When selecting Compound as the Network Type, a dropdown menu appears under Query, providing the option of choosing between compounds or genes. When any other Network Type is selected, the only option is to use compounds/genes unless using a selected pathway.

10. Click **Build Network** to query the database and create the network.

Option 2: Load an experiment file

Use to load experimental data to visualize and explore compound networks over time or in varying experimental conditions. The input can be an Excel, comma, or tab delimited file. MetScape allows users to load three types of files – compound file, gene file, and concept file. Each type is optional, e.g. you can load only compounds, only genes, only concepts, or any combination of the above.

Compound File

The compound data file must meet the following requirements:

- The first row must be a heading row, in which:
 - The columns in the first row are column headings to label the data.
- All other rows contain experimental data, in which:
 - The first column contains KEGG Compound IDs or names.
 - The remaining columns contain experimental data. Multiple experimental values are permissible in the same spreadsheet.

Below is a portion of an example compound file with Compound IDs, significance values, and fold change values:

	A	B	C
1		Fold Change	P-value
2	C00002	0.51	0.006174265
3	C00004	1.14	0.646801757
4	C00008	0.99	0.970303452
5	C00020	2.12	0.068052875
6	C00025	0.91	0.115376267
7	C00037	0.99	0.904658012
8	C00041	0.95	0.574835153
9	C00042	1.09	0.809172269
10	C00044	0.56	0.100724407
11	C00047	1.36	0.189989302
12	C00049	0.92	0.755789133
13	C00064	0.48	0.186707182
14	C00065	1.22	0.805151078
15	C00073	0.80	0.356405722

Gene File

The gene file must meet the following requirements:

- The first row must be a heading row that includes column headings to label the data.
- The first column contains Entrez Gene IDs or Official Gene Symbols.
- The remaining columns contain experimental data. Multiple experimental values are permissible in the same spreadsheet.

Below is a portion of an example gene file with Gene IDs, significance values, and fold change values:

	A	B	C
1	ENTREZ_GENE_ID	Pvalue	Log_fold_change
2	1	0.002335904	-0.02
3	2	2.07E-08	-0.06
4	2	2.03E-08	0.04
5	9	0.003919472	0.04
6	10	0.00026716	-0.08
7	12	0.000237925	0.05
8	13	0.002489364	-0.05
9	14	0.506396205	0.00
10	15	8.25E-06	-0.04
11	16	0.022222108	-0.01
12	18	1.47E-08	-0.07
13	18	5.24E-06	-0.09
14	18	2.12E-05	-0.07
15	19	0.820594937	0.00
16	19	1.76E-08	0.06
17	19	9.46E-09	0.04
18	19	0.660631201	0.00
19	20	1.16E-07	-0.05


Concept File

The concept file can be generated by a gene set enrichment analysis tool such as LRpath or GSEA from gene expression data.


need a Gene Expression file with fold change (or log fold change) values and p-values. The Gene Expression file needs to contain all gene records, not just those for significant genes; LRpath will determine the significant genes from the input.

5. Click **OK**.

Example of selecting experimental data, including importing files and designating fold change, p-values, and thresholds:




Example of ID mapping window from a loaded compound experiment file:



Note: If you select **Save** on the **Select Compound Mappings** window, you will get a .csv file containing mapping information for all compounds that MetScape successfully mapped.

Note: The **Select Experimental Data** window can also be accessed from the MetScape tab.

Note: Multiple networks can be built within the same user session. After a network is built, pull up the **Select Experimental Data** window and change the data files. Then build the new network. Both networks will remain accessible. If you change species, however, all existing data will be lost. The below warning window will pop up in this situation.



Missing Data Window


After the experimental data has been loaded, if any of the genes/compounds that you submitted were not mapped to the database objects, a **Missing Data** window will appear.

Genes, compounds, and concepts may appear on this missing elements list because:

- Genes and compounds that you supply may not be in the database. If they are not found in the database then they are reported as missing.


- If your input genes are not human (Rat, for example), then they are mapped to human genes using homologs from NCBI's HomoloGene. If this mapping fails, then those genes are reported as missing.
- MetScape will display only the genes that encode metabolic enzymes. If an input gene does not encode metabolic enzymes, it will appear on the missing elements list.
- A concept (pathway) will appear on the missing list if all of its significant genes are missing. The list of significant genes for a concept comes from the input file or from LRpath.

After viewing the **Missing Data** information, click **SAVE** to save the data for later viewing or click **OK** to close the window without making it available for future viewing.



VISUALIZATION

Access the MetScape **Legend** from **Plugins -> MetScape -> Show Legend**.



Additional information about the network is expressed through visualization:

- A compound is red in the original network and subsequent subnetworks if it was in the original data loaded into MetScape.
- A green border surrounding a node represents significant gene/compound.
- Node size represents the direction of the change. Larger nodes represent an increase and smaller nodes represent a decrease in gene/metabolite. The actual amount of the change is not represented visually.
- When a node is expanded, the edges between the original node and the expanded nodes become blue.

METSCAPE TAB

This tab has options for choosing Network Types, the data you entered, and more.

From this tab you can:


- **Select the type of network that you want to build**

Select **Network Type** by selecting one of the following from the dropdown menu:


- **Compound-Reaction-Enzyme-Gene**
- **Compound-Reaction**
- **Compound-Gene**
- **Compound**

Note: When selecting Compound as the Network Type, a dropdown menu appears under Query, providing the option of choosing between compounds or genes. When any other Network Type is selected, the only option is to use compounds/genes unless using a selected pathway.

- **Add or Remove** data.
- Select a **pathway**
 - Use to view all the compounds and reactions associated with a metabolic pathway. Select a specific pathway from the drop-down list.



- Save the **Output as a File**.
- **Build Network** graph based on the data.



Note: When the organism is human, the Genes section will show Input ID and Input Symbol. However, when the organism is not human (ex. rat), the Genes section will show Input ID, Input Symbol, and Human Symbol.

RULES USED TO BUILD DIFFERENT NETWORK TYPES

- **Compound-Reaction-Enzyme-Gene (C-R-E-G), Compound-Reaction (C-R), and Compound-Gene (C-G) networks**

The C-R-E-G, C-R, and C-G networks are all built from the same underlying data. That data is derived, in each case, by finding compounds that participate in reactions that are catalyzed by enzymes that are encoded by genes. If only genes are input, then all the enzymes, reactions, and compounds that match those genes are used. If only compounds are input, then all the reactions, enzymes, and genes that match those compounds are used. If both genes and compounds are input, then only those C-R-E-G couplings that match both a compound from the compound input and a gene from the gene input are used.

If a concept file is provided, genes from that file will be used as input. If a gene file is also provided it will be used only as input for LRpath (if user chooses that option).

If only a gene file is provided, all genes from that file are used as input. In this case, we recommend that you load a smaller set of genes (e.g. most significant differentially expressed genes).

- **Compound (C) networks built from compound input**

If a list of compounds is provided, the resulting network will include the query compounds (shown in red) plus any compounds that participate in the same reactions as query compounds. The edges will be drawn between “seed” compounds and their neighboring compounds

- **Compound (C) networks built from gene input**

If a list of genes is provided, the resulting network will include the compounds that are related to the query genes via the reactions in which they participate, and the enzymes that catalyze these reactions. The edges will be drawn between all compounds.

- **Pathway networks**

Pathway-specific C-R-E-G, C-R, and C-G networks are built from a set of genes, enzymes, reactions and compounds defined in the EHMN database. Pathway-specific Compound networks are built similarly to the Compound networks built from a set of input genes.


DATA PANEL DISPLAY

Attributes of compounds, reactions, pathways, and concepts that you have selected in the network appear in the Cytoscape Data Panel.


Choosing attributes provided through MetScape

To choose which attributes to view:

1. To choose attributes for compounds, click the **Node Attribute Browser** tab in the Data Panel. You can choose to view **node**, **edge**, or **network** attributes by clicking on the respective tabs at the bottom of the Data Panel. Pathway data can be viewed under the **Pathway Filter tab** and concept data can be viewed under the **Concept Filter** tab in the Data Panel.
2. Click the **Select Attributes** icon in the Data Panel toolbar. A list of attributes with check boxes will appear.



3. Select the attributes you want displayed in the Data Panel. Attributes are listed with item type before the attribute (ex. Compound.name; Enzyme.name). The available attributes will depend on the type of network created.



4. When you are done, click anywhere outside of the list and the table will be populated with the appropriate data.

Below is a table of available attributes for each node type:

A	B
Attribute	Description
ID	Unique identifier, required by Cytoscape
Compound.approximateMW	Compound Molecular weight (integer, data source - EHMN database)
Compound.casnum	Chemical Abstracts Service compound identifier
Compound.cid	KEGG compound unique identifier
Compound.file_name	Rows with data specific to your experiments
Compound.formula	Molecular formula
Compound.mw	Molecular weight
Compound.name	Compound name
Enzyme.ecnum	Enzyme Comission number
Enzyme.name	Enzyme name
Gene.file_name	Relative gene data values
Gene.column_name	Gene data p-value
Gene.description	Gene description
Gene.human.geneid	Human gene identifier
Gene.symbol	Official gene symbol
Gene...geneid	Ortholog gene ID if the input organism is not human
Reaction.equation	Equation for a reaction
Reaction.locations	Reaction subcellular location (EHMN database)
Reaction.pathway	Pathway involving reaction
Reaction.reversible	T= reaction is reversible; F=reaction is non reversible
Reaction.rid	Reaction identifier
Type	Node type (gene, enzyme, reaction or compound)
canonicalName	Node name that is shown on node label
Gene.locations	Subcellular location

*Please note that the reactions shown in the Reaction.equation attribute display all compounds, while the networks built from the list of input genes/compounds show only "main" compounds (see [Ma et al., Mol. Syst. Biol., 2007, 3:135](#) for details). Pathway-specific networks show all compounds.

Below is a table of available attributes for each edge:

A	B
Attribute	Description
Enzyme.names	Enzyme name
Enzyme.ecnums	Enzyme Comission number
Reaction.equation	Equation for reaction
Reaction.locations	Reaction subcellular location (EHMN database)
Reaction.pathway	Pathway involving reaction
Reaction.reversible	T= reaction is reversible; F=reaction is non reversible
Reaction.rid	Unique identifier for reaction
canonicalName	Common name
interaction	Type of interaction

5. Select compounds and/or edges in the graph to view their attributes in the Data Panel.

ID	Compound.formula	Compound.name
R01218		
CE4852	C43H64N7O17P3S	trans,cis,cis,cis,cis-2,10,13,16,19-docosapentaeno...
516		
1.1.99.2		
C05108	C29H48O	14-Demethylsterol
R01664		
RE3553		
131		

Node Attribute Browser | Edge Attribute Browser | Network Attribute Browser | Concept Selection

Click on heading and drag with mouse (notice: Compound.name is currently to the right of Compound.formula)

ID	Compound.name	Compound.formula
R01218		
CE4852	trans,cis,cis,cis,cis-2,10,13,16,19-docosapentaeno...	C43H64N7O17P3S
516		
1.1.99.2		
C05108	14-Demethylsterol	C29H48O
R01664		
RE3553		
131		

Node Attribute Browser | Edge Attribute Browser | Network Attribute Browser | Concept Selection

Drop heading in new location (notice: Compound.name is now next to ID and to the left of Compound.formula)

To sort by an attribute, click the attribute column heading; to sort in the reverse direction, click again.


ID	Compound.name ▾	Compound.formula
CE4814	(8Z,11Z,14Z,17Z)-cisoctadecenoic OJA	C41H62N7O17P3S
G00119	(Gal)1 (GalNAc)1 (Glc)1 (Neu5Ac)3 (Cer)1	
G00120	(Gal)2 (GalNAc)1 (Glc)1 (Neu5Ac)3 (Cer)1	
C01205	(R)-3-Amino-2-methylpropanoate	C4H9NO2
C04677	1-(5'-Phosphoribosyl)-5-amino-4-imidazolecarboxa...	C9H15N4O8P
C04734	1-(5'-Phosphoribosyl)-5-formamido-4-imidazolecar.	
CE5944	10,11-dihydro-12-oxo-LTB4	
C00234	10-Formyltetrahydrofolate	

Node Attribute Browser | Edge Attribute Browser | Network Attribute Browser

Click on column title to sort. Notice the ▾ (arrow) appears when a column is sorted A-Z (an arrow in the reverse direction appears when sorted Z-A).

Pathway filter in the data panel

MetScape provides two ways to access pathway information. The **Pathway Filter** tab lists all pathways represented in the network. Selecting one or more pathways will highlight in the network all nodes in the pathway(s). Pathways are also displayed as attributes in the data panel.



A subgraph can be created for a given pathway. The subgraph consists of a subset of nodes from an active network that belong to the selected pathway. To create a subgraph:

1. Select a pathway in the **Pathway Filter** tab.
2. Click the **Create Subnetwork** button at the top of the Data Panel (make sure you are on the **Pathway Filter** tab).
3. A new graph appears in the graph window. This is the subnetwork for the selected pathway.

Concept filter in the data panel

The content of the created or loaded concept file is put into the **Concept Filter** tab.

The screenshot shows the 'Data Panel' window with the 'Concept Filter' tab selected and highlighted with a red box. The table below displays enrichment data for various concepts. Arrows from external labels point to specific columns: 'Concept names' points to the 'Concept Name' column, 'Direction of the change' points to the '↑↓' column, 'P-value' points to the 'P-value' column, and 'FDR' points to the 'False Discovery Rate' column.

Concept Name	Number of Enriched-Driving Genes	Number of Mapped Genes	↑↓	P-value	False Discovery Rate
Alzheimer's disease	127	9	↑	2.56E-6	3.64E-5
Base excision repair	26	3	↑	3.08E-4	0.00229
Biosynthesis of alkaloids derived from histidine and purine	51	12	↑	0.00108	0.00636
Biosynthesis of alkaloids derived from ornithine, lysine and nicotinic acid	45	7	↑	0.00686	0.0293
Biosynthesis of alkaloids derived from shikimate pathway	44	7	↑	0.0137	0.0497
Biosynthesis of alkaloids derived from terpenoid and polyketide	55	17	↑	6.31E-11	3.6E-9
Biosynthesis of phenylpropanoids	41	5	↑	0.00408	0.0194
Biosynthesis of plant hormones	64	18	↑	7.53E-9	1.84E-7
Biosynthesis of terpenoids and steroids	50	15	↑	5.72E-9	1.63E-7

1. Concept Name - official name of the concept.
2. Number of Genes - the number of input genes that belong to the concept.
3. Direction - direction of the change.
4. P-value - the p-value for enrichment/depletion.

5. False Discovery Rate - the FDR, estimated using the method of Benjamini and Hochberg (When results are sorted by p-value, a FDR < 0.05 is interpreted to mean that approximately 5% of the concepts with lesser or equal p-value are false positives).

To map a concept on the graph, click on the concept name in the Data Panel. The concept will then be highlighted in the graph (associated nodes turn yellow and associated edges turn red).

Name	Number of Enriched-Driving Genes	Number of Mapped Genes	↑↓	P-value	False Discovery Rate
Biosynthesis of plant hormones	64	18	↑	7.53E-9	1.84E-7
Biosynthesis of terpenoids and steroids	50	15	↑	5.72E-9	1.63E-7
Biosynthesis of saturated fatty acids	19	4	↑	0.00867	0.0337
Cell cycle	100	2	↑	1.6E-5	1.95E-4
Citrate cycle (TCA cycle)	25	8	↑	4.78E-5	4.3E-4
Cytokine-cytokine receptor interaction	132	6	↓	3.05E-7	5.8E-6
DNA replication	30	3	↑	5.46E-11	3.6E-9
Focal adhesion	153	4	↓	1.44E-4	0.00117
Glycine, serine and threonine metabolism	23	2	↓	0.00728	0.0296

A subgraph can be created for a given concept. The subgraph consists of genes considered to be significant to the concept (those genes that drive the concept), not all genes in the concept. To create a subgraph:

4. Select a concept in the **Concept Filter** tab.
5. Click the **Create Subnetwork** button at the top of the Data Panel (make sure you are on the **Concept Filter** tab).
6. A new graph appears in the graph window. This is the subnetwork for the selected concept.

Graph of subnetwork appears in graph window

Select concept


Click on "Create Subnetwork" button

Concept Name	Number of Enriching Genes	Number of Mapped Genes	↑↓	P-value	False Discovery Rate
Biosynthesis of plant hormones	64	18	↑	7.53E-9	1.84E-7
Biosynthesis of terpenoids and steroids	50	15	↑	5.72E-9	1.63E-7
Biosynthesis of unsaturated fatty acids	19	4	↑	0.00867	0.0337
Cell cycle	100		↑	1.6E-5	1.95E-4
Citrate cycle (TCA cycle)	25		↑	4.78E-5	4.3E-4
Cytokine-cytokine receptor interaction	132		↑	3.05E-7	5.8E-6
DNA replication	30		↑	5.46E-11	3.6E-9

The **Reapply Selection** button will re-select the last selected concept in the graph. If you select a concept and then select off of it, clicking the **Reapply Selection** button will re-select that concept.

The concept information can be saved as a file. To save this information:


1. Under **Concept Filter** tab, click the **Save Concepts...** button.
2. A new window pops up, asking where to save the file. Choose the desired location.
3. Click **Save**. The file saves as a CSV file that can be opened with most text editing or spreadsheet programs. The file format is identical to the input concept file.




METSCAPE RESULTS PANEL (ADDITIONAL NODE INFORMATION)

Additional information can be obtained for each node by double clicking on it.

To get additional information about any compound, gene, reaction, or enzyme:

1. Double click on the node of interest.
2. On the right side of the screen, the **MetScape Results Panel** will show up.
3. Use the scroll bar or click on the undock icon  to see all the data.



ANIMATION OF DATA


Create an animation of the data to see how it changes over time and across treatments.

Building an Animation

1. Create a subnetwork (for example, create a subnetwork for the TCA Cycle. See the [Concept Selection Data](#) in the Data Panel section for how to do this).


Note: It is not necessary to create a subnetwork before doing data animation; the animation can be done on any network.

2. From the **Plugins** menu, choose **MetScape -> Animate data...**




3. A new window pops up, titled **Initialize Coupled Animation of Multiple Data Columns**.
4. Column headings from your data are listed on the right side of the new window. Determine how many animations (number of rows) and data items per animation

(number of columns) you will need. Enter this information into the Rows and Cols boxes. If you change these numbers from the default, click on **Reconfigure Layout**.




5. Click and drag the column names into the grid to build the animation.



6. After grid is completely filled in, click on **Build Animation**.

Note: If you want to undo what you filled out in the grid, click on **Restart** and the grid will be cleared.

7. Three new windows appear in the graph panel, each representing a different treatment. In addition, an **Animation Controls** window appears.



Manipulating the animation: zooming

1. **Zoom** in on one graph to the appropriate view size (zoom the same way you do with any Cytoscape graph).
2. On the **Animation Controls** window, click on **Realign All**. This brings all the treatment graphs to the same view.


4.2.1.3
3420
C00311
50


Zoom in on one graph

Click Realign All button

The other 2 graphs align to match the zoomed in graph

Manipulating the animation: bar chart

1. The bar chart in the **Animation Controls** window shows the range of the data and its frequency. The range of measurements are on the x-axis, while the frequency with which they occur in the experimental data set are on the y-axis.
- 
2. Click on a bar in the bar chart and a vertical line (slider) will show up on the bar chart. The slider can be moved by clicking on the bar chart and dragging the line. When moved, the values shown in black are also shown in black in the graph and the set of the color scale is reset.
 3. The selected portion of the graph determines what the animation will show. This allows outliers to be removed if desired.
 4. The color range can be changed by clicking the dropdown arrow next to **Red-Blue**.



Playing the animation

To play the animation, click **Play** on the **Animation Controls** window. Colors change within the three treatments, showing changes over time. To stop the animation, click **Stop** on the **Animation Controls** window.

The screenshot displays the MetScape software interface. At the top, three network windows are shown: DataSeries_01, DataSeries_11, and DataSeries_21. Each window contains a metabolic network diagram with nodes and edges. The nodes are color-coded: blue, green, orange, and purple. The Animation Controls window is visible at the bottom, featuring a 'Stop' button and a 'Close' button. A yellow callout box with black text reads 'Colors change representing change over time', with arrows pointing to the nodes in the three network windows.


MANIPULATING NETWORKS

Expanding a network

From a compound node, you can expand the current network to include additional reactions and related compounds, genes, and enzymes.


Option 1: Expanding a network in the current network window (Example using Compound-Reaction Network)

1. Right click on the compound node you want to expand (this node is known as the expansion seed node). A menu of options will pop up.
2. Go to **MetScape** -> **Expand** -> **Expand in Existing Network**




3. Additional compounds and reactions are added to the network. As a result, the network is often redrawn.

4. The edges between the original node and the expanded nodes are now blue.




Option 2: Expanding a network in a new window (Example using Compound Network from Selected Genes)

1. Right click on the compound node you want to expand (this node is known as the expansion seed node). A menu of options will pop up.
2. Select **MetScope -> Expand -> Expand in Subnetwork**



3. A new network is created in a new window. This network includes only the expansion seed compound and its related compounds, reactions, enzymes, and genes (depending on the Network Type). The original color designations are used when the graph is created in a new window.



Collapsing a network

To collapse a network that is expanded in a current network window:


1. Right click on the expansion seed compound node.
2. Select **MetScape -> Collapse** and then:
 - To collapse only the branch expanded from that expansion seed compound, select **Collapse**.
 - To collapse all expanded branches, select **Restore Original Network**.

Creating a subnetwork

A subnetwork of a current network can be created and will appear in a separate window. A subnetwork will include all highlighted nodes and edges (selected nodes should be yellow and selected edges should be red).

For example:


1. Select the compound IMP
2. Go to the **Select** menu in Cytoscape.
3. Choose **Nodes -> First Neighbors of Selected Nodes**. Now all the first neighbors of IMP should be yellow.
4. Go to the **Select** menu in Cytoscape.
5. Choose **Edges -> Select Adjacent Edges**. Now all the adjacent edges should be red.
6. Right click on IMP. A menu of options will pop up.



7. Choose **MetScape -> Create Subnetwork**.

Note: When choosing options from the menu, be careful not to move the cursor outside the menu panels. Doing so will remove the highlighting of the nodes and edges, resulting in an empty subnetwork.

8. A new subnetwork is created in the graph window.



Destroying a network

To destroy a network no longer needed:

1. Make sure you really do want to destroy the network.
 - Cytoscape will not ask you to confirm the deletion.
 - Destroying the network is irreversible.
 - **Note:** If you want to destroy a network view without destroying the network itself, use **Destroy View** instead of **Destroy Network**.
2. On the **Network** tab in the Cytoscape Control Panel, right click on the network you want to destroy and select **Destroy Network**.

Note: If you destroy a network that has subnetworks, the subnetworks are **NOT** destroyed. They are promoted up one network level.

SAVING AND REOPENING A SESSION

Saving a session

To save a Cytoscape session containing one or more MetScape plugin networks:

1. Select **File -> Save** (or **Save As...**) from the Cytoscape menu.
2. Browse to a location for saving the file.
3. Name the file.
4. Click **Save**.

Reopening a session

To reopen a saved session containing MetScape plugin data:

1. Select **File -> Open** from the Cytoscape menu.
2. Navigate to the saved file location.
3. Select the file.
4. Click **Open**.
5. After Cytoscape reports that the session file was successfully loaded, click **Close**.